

Aurora Discovery TCPR



Manufacturer: Aurora Discovery
Model Number: tcPR
Web Address: www.auroradiscovery.com

The Topology Compensating Plate Reader (tcPR) is the second generation emission-ratioing, fluorescence reader designed to read 1536-well and 3456-well assay plates. It acquires fluorescence from each well of a plate simultaneously at two different wavelengths, thus providing ratioing capability. Excitation is in the range 340-650 nm using a xenon lamp light source. Emitted light is split between two optical fibers, passed through two emission filters, each at different wavelengths (in the range of 400-600 nm), and collected by photomultiplier tubes.

Topology Compensating Plate Reader: Commands

- ▶ **PlateIn()** - Command the tcRP to lower (take in) its plate carrier.
- ▶ **PlateOut()** - Command the tcRP to present its plate carrier for loading.
- ▶ **SetScriptFile(path)** - Set the path name for the script file.

path	String	Full path name for the script file.
------	--------	-------------------------------------

- ▶ **SetPlateType(plate)** - Set the name of the plate type to be used for the read.

plate	String	Name of the plate type.
-------	--------	-------------------------

- ▶ **ReadPlate()** - Run the read operation.
- ▶ **GetVersionString()** - Returns the string containing the version number of the instrument.

Topology Compensating Plate Reader: Errors

- ▶ **ReaderError(errorCode, errorDescription)** - Error occurred during command execution.

errorCode	Integer	Error code number.
-----------	---------	--------------------

errorDescription	String	Error description. 1 - Toolkit is already initialized. 2 - Product path specified in the Initialize method is invalid. 3 - Toolkit fails to initialize the tcPR. 4 - Bad command. 5 - Bad parameter. other - Please refer to tcPR error data file.
------------------	--------	--

ReTiSoft Inc.
366 Revus Avenue, Unit 21
Mississauga, Ontario, Canada, L5G-4S5
Main: 647-724-2398 Europe: 33-9-7518-0225
Web: www.retisoft.ca Email: prodziew@retisoft.ca