

QInstruments BioShake 3000 T-ELM



Manufacturer: QUANTIFOIL Instruments GmbH
 Model Number: QInstruments-Bioshake3000T-ELM
 Web Address: www.qinstruments.com

Bioshake3000-T-ELM is the world's smallest industrial All-in-One ThermoShaker. It is specially designed for integration in robotic environments. It provides Automatic Edge Locking Mechanism (ELM) for robotic gripping, fast and high-precision heating from ambient to 99C, mixing speed up to 3,000 rpm at 2 mm orbit, and homing position accuracy of +/- 0.1mm.

Shaker: Commands

- ▶ **resetDevice()** - Restart the controller.
- ▶ **shakeOn()** - Start the shaking with the current mixing speed.
- ▶ **shakeOff()** - Stop the shaking and return to homing position.
- ▶ **getShakeState()** - Return the state of the shaker.
- ▶ **getShakeTargetSpeed()** - Return the target mixing speed.
- ▶ **setShakeTargetSpeed(speed)** - Set the target mixing speed.

| | | |
|-------|---------|----------------------|
| speed | Integer | Target mixing speed. |
|-------|---------|----------------------|

- ▶ **getShakeActualSpeed()** - Return the current mixing speed.
- ▶ **getShakeMinRpm()** - Return the least set point.
- ▶ **getShakeMaxRpm()** - Return the biggest set point.
- ▶ **tempOn()** - Activate the temperature control.
- ▶ **tempOff()** - Deactivate the temperature control.
- ▶ **getTargetTemp()** - Return the target temperature.
- ▶ **setTargetTemp(temperature)** - Set the target temperature in 1/10 degrees.

| | | |
|-------------|---------|---------------------|
| temperature | Integer | Target temperature. |
|-------------|---------|---------------------|

- ▶ **getActualTemp()** - Return the current temperature.
- ▶ **getTempMin()** - Return the least set point of temperature.
- ▶ **getTempMax()** - Return the biggest set point of temperature.
- ▶ **getTemp40Calibr()** - Return the temperature on the calibration point at 40 degree Celsius.
- ▶ **setTemp40Calibr(temperature)** - Set the temperature on the calibration point at 40 degree Celsius in 1/10 degrees.

| | | |
|-------------|---------|-----------------------------------|
| temperature | Integer | Temperature on calibration point. |
|-------------|---------|-----------------------------------|

▶ **getTemp90Calibr()** - Return the temperature on the calibration point at 90 degree Celsius.

▶ **setTemp90Calibr(temperature)** - Set the temperature on the calibration point at 90 degree Celsius in 1/10 degrees.

| | | |
|-------------|---------|-----------------------------------|
| temperature | Integer | Temperature on calibration point. |
|-------------|---------|-----------------------------------|

▶ **setTempLoggingOn()** - Enable the temperature logging function.

▶ **setTempLoggingOff()** - Disable the temperature logging function.

▶ **setElmShakePos()** - Close the Edge Locking Mechanism (ELM). The microplate is now fixed.

▶ **setElmUnlockPos()** - Opens the Edge Locking Mechanism (ELM) for gripping microplates.

▶ **getElmState()** - Return the state of Edge Locking Mechanism (ELM).

BioShake Error

🔥 **BioShakeError(code, description)** - Error occurred during command execution.

| | | |
|------|---------|-------------|
| code | Integer | Error code. |
|------|---------|-------------|

| | | |
|-------------|--------|--------------------|
| description | String | Error description. |
|-------------|--------|--------------------|

ReTiSoft Inc.
 366 Revus Avenue, Unit 21
 Mississauga, Ontario, Canada, L5G-4S5
 Main: 647-724-2398 Europe: 33-9-7518-0225
 Web: www.retisoft.ca Email: info@retisoft.ca