



# Präzisionstechnik

SPANNTECHNIK / AUTOMATION

ROBOTIK / FÖRDERTECHNIK

## Programmable Servo Gripper for the Laboratory Automation



The PTM Servo Gripper with serial interface optimizes the sensitive handling of the various products and devices in the laboratory automation as well as applications in the general automation industry area

- Pipettes
- Cuvettes
- Tip Boxes
- Test Tubes
- Microplates
- Microplate-Lids

P.T.M. Präzisionstechnik GmbH  
Olchinger Straße 109  
D 82194 Gröbenzell  
Germany  
Tel. +49 (0)8142 59 39 7-0  
Fax +49 (0)8142 59 39 7-28  
E-mail: info@ptm-automation.de

The wide stroke of the gripper provides handling of microplates in both orientations without change of fingers

