


Beckman Biomek 2000



Manufacturer: Beckman Coulter, Inc.
Model Number: Biomek2000
Web Address: www.beckmancoulter.com

The Biomek 2000 Laboratory Automation Workstation meets the demands of rapidly changing life science technology with simple, intelligent automation of liquid-handling tasks. Pipetting, diluting and dispensing operations are performed quickly, easily and automatically. The modular platform allows expansion of system capability to include plate heating and cooling, plate washing, high-density transfers, photometric measurement and high-capacity operation. The entire system is controlled by powerful and user-friendly BioWorks software with an intuitive graphical interface.

Biomek 2000 Workstation: Commands

▶ **Execute(method)** - Start executing the specified method on the Biomek 2000.

| | | |
|--------|--------|---|
| method | String | Name of the method to execute on the Biomek 2000. |
|--------|--------|---|

Biomek 2000 Workstation: Errors

🔥 **BiomekError(code)** - Error occurred during command execution.

| | | |
|------|--------|-------------|
| code | String | Error code. |
|------|--------|-------------|

ReTiSoft Inc.
 366 Revus Avenue, Unit 21
 Mississauga, Ontario, Canada, L5G-4S5
 Main: 647-724-2398 Europe: 33-9-7518-0225
 Web: www.retisoft.ca Email: prodziew@retisoft.ca