

Beckman Multimek



Manufacturer: Beckman Coulter
 Model Number: Beckman-Multimek
 Web Address: www.beckmancoulter.com

The Multimek 96/384-Channel Automated Pipettor from Beckman Coulter has the flexibility to transfer between 96-, 384- and 1536-well plates on the same worksurface. The Multimek 96/384 delivers high performance and reproducibility while meeting the applications requirements of today's laboratory.

Multimek: Commands

▶ **LoadScript(script)** - Load a script into memory.

script	String	Specifies the script file to load.
--------	--------	------------------------------------

▶ **RunScript()** - Run the currently loaded script.

▶ **GetStatus()** - Get the status of the instrument, one of ERROR, PAUSED, READY, RUNNING, or WAITING.

▶ **GetError()** - Get the cause of the error.

▶ **ResetError()** - Reset the error.

Multimek: Errors

🔥 **PipettorError(code, alarmMessage)** - Error occurred during command execution.

code	Integer	Error code.
alarmMessage	String	Error message for the instrument

ReTiSoft Inc.
 366 Revus Avenue, Unit 21
 Mississauga, Ontario, Canada, L5G-4S5
 Main: 647-724-2398 Europe: 33-9-7518-0225
 Web: www.retisoft.ca Email: info@retisoft.ca