

(c) ReTiSoft Inc.

www.retisoft.ca

Cognex DM100QL



Manufacturer: Cognex Systems, Inc.
Model Number: Cognex-DM100QL
Web Address: www.cognex.com

The DataMan 100 Series readers offer industry-leading reading performance in a small product size with USB and RS-232 communication options.

Scanner: Commands

- **TriggerOn()** - Turns the barcode laser on and off with a timeout.
- **TriggerOff()** - Turns the barcode laser off.

Barcode-Reader Error

🚫 **BarcodeError(alarmCode, description)** - Error occurred during command execution.

alarmCode	Integer	Error code indicating an error which occurred on the instrument.
description	String	Error description.

ReTiSoft Inc.
366 Revus Avenue, Unit 21
Mississauga, Ontario, Canada, L5G-4S5
Main: 647-724-2398 Europe: 33-9-7518-0225
Web: www.retisoft.ca Email: info@retisoft.ca