

CyBio CyBi-Well



Manufacturer: CyBio AG
Model Number: CyBi-Well
Web Address: www.cybio-ag.com

CyBi-Well can be optionally configured with a 96-well or 384-well pipetting head. With its 96 or 384 simultaneously operating precision pipettes, the pipettor is ideal for large or complex test series in the HTS and uHTS range. In only 20 seconds it pipettes and transports a 384-well microplate; it requires 30 seconds and four steps to process 1,536-well microplates. The highly precise lifters enable the use of shallow-well microtiter plates in deep-well devices.

Liquid Handler: Commands

▶ **OpenDocument(path)** - Opens the specified document.

path	String	Full file-path of the document to be opened.
------	--------	--

▶ **ExecuteDocument(path)** - Executes the specified document.

path	String	Full file-path of the document to be executed.
------	--------	--

▶ **CloseDocument(path)** - Closes the specified document.

path	String	Full file-path of the document to be closed.
------	--------	--

Liquid Handler: Errors

🔥 **LiquidHandlerError(code, message)** - Error occurred during command execution.

code	Integer	Error code.
------	---------	-------------

message	String	Error message.
---------	--------	----------------

ReTiSoft Inc.
 366 Revus Avenue, Unit 21
 Mississauga, Ontario, Canada, L5G-4S5
 Main: 647-724-2398 Europe: 33-9-7518-0225
 Web: www.retisoft.ca Email: prodziew@retisoft.ca