

MediaCybernetics Image-Pro



Manufacturer: MediaCybernetics
Model Number: Auto-Pro
Web Address: www.mediacy.com

Image-Pro Plus is the ultimate image analysis software package for fluorescence imaging, quality assurance, materials imaging, and various other scientific, medical, and industrial applications. The most powerful member of the Image-Pro software family, Image-Pro Plus includes extensive enhancement and measurement tools and allows you to write application-specific macros and plug-ins.

Auto-Pro Software: Commands

▶ **Scan(barcode)** - Executes the Auto-Pro functions to scan the plate.

barcode	String	barcode of the plate.
---------	--------	-----------------------

Auto-Pro Software: Errors

🔥 **Error(number)** - Error occurred during command execution.

number	Integer	Error number.
--------	---------	---------------

ReTiSoft Inc.
 366 Revus Avenue, Unit 21
 Mississauga, Ontario, Canada, L5G-4S5
 Main: 647-724-2398 Europe: 33-9-7518-0225
 Web: www.retisoft.ca Email: prodziew@retisoft.ca