

## Microscan VS310



**Manufacturer:** Microscan Systems, Inc.  
**Model Number:** VS310  
**Web Address:** [www.microscan.com](http://www.microscan.com)

Scanner capable of scanning Code 39, Codabar, Interleaved 2 of 5, Code 128 and Code 93 types of bar codes.

### Scanner: Commands

▶ **r( )** - Turns the barcode laser on. If the barcode can not be read correctly the "NOREAD" string is returned.

ReTiSoft Inc.  
366 Revus Avenue, Unit 21  
Mississauga, Ontario, Canada, L5G-4S5  
Main: 647-724-2398 Europe: 33-9-7518-0225  
Web: [www.retisoft.ca](http://www.retisoft.ca) Email: [prodziew@retisoft.ca](mailto:prodziew@retisoft.ca)