

PerkinElmer Evolution P3



Manufacturer: PerkinElmer Life Sciences
Model Number: Evolution P3
Web Address: www.perkinelmer.com

EvolutionP3 offers a 96 or 384 channel Dispense Head for dispensing into 96, 384 and high density 1536 well plates, a 16 position labware deck and an automatic tip load feature. A wide dynamic volume range enables high precision pipetting from 0.5ul, and up with disposable tips and dispense volumes as low as 20ul using Pintool technology.

Liquid Handler: Commands

▶ **execute(name)** - Runs a script to operate the P3 liquid handler.

name	String	Specifies the script name and its location.
------	--------	---

ReTiSoft Inc.
 366 Revus Avenue, Unit 21
 Mississauga, Ontario, Canada, L5G-4S5
 Main: 647-724-2398 Europe: 33-9-7518-0225
 Web: www.retisoft.ca Email: prodziew@retisoft.ca