

PerkinElmer JANUS



Manufacturer: PerkinElmer Life Sciences
Model Number: JANUS
Web Address: www.perkinelmer.com

The JANUS automated workstation is designed for the efficient automation of sample preparation procedures utilized in pharmaceutical, biotech and research applications. The instrument has a modular platform which consists of one pipetting arm with different tip configurations as well as one plate movement arm.

Liquid Handler: Commands

▶ **execute(name)** - Runs a script to operate the JANUS liquid handler.

name	String	Specifies the script name and its location.
------	--------	---

ReTiSoft Inc.
366 Revus Avenue, Unit 21
Mississauga, Ontario, Canada, L5G-4S5
Main: 647-724-2398 Europe: 33-9-7518-0225
Web: www.retisoft.ca Email: prodziew@retisoft.ca