

PerkinElmer MultiPROBE II EX (Left and Right) with Gripper



Manufacturer: PerkinElmer Life Sciences
Model Number: MultiPROBE II EX
Web Address: www.perkinelmer.com

The Gripper Integration Platform is a powerful device that expands the power of 4-tip and 8-tip MultiPROBE II EX systems. The Gripper Platform provides a secondary arm with integrated gripper tool capable of picking-and-placing SBS-approved microplates, microplate lids, deep-well plates, extraction blocks and selected vacuum manifolds around the deck of MultiPROBE II EX systems.

Liquid Handler: Commands

▶ **execute(name)** - Runs a script to operate the MultiPROBE II EX liquid handler.

name	String	Specifies the script name and its location.
------	--------	---

ReTiSoft Inc.
 366 Revus Avenue, Unit 21
 Mississauga, Ontario, Canada, L5G-4S5
 Main: 647-724-2398 Europe: 33-9-7518-0225
 Web: www.retisoft.ca Email: prodziew@retisoft.ca