

PerkinElmer Victor



Manufacturer: PerkinElmer Life Sciences
Model Number: Victor
Web Address: www.perkinelmer.com

Victor multilable counter is a complete platform for quantitative detection of light-emitting, or light-absorbing markers. It is suitable for flush or glow luminometry, UV absorbance, photometry, homogeneous time-resolved fluorescence assays and fluorescence polarization. It is a compact bench top unit with features such as dispensers, shaker, temperature control, top or bottom counting or scanning.

Plate Reader: Commands

▶ **RunAssay(protocolNumber)** - Load the plate into the reader and execute the assay.

| | | |
|----------------|---------|--|
| protocolNumber | Integer | protocol number for the assay. Range of values: 1 to 99 . |
|----------------|---------|--|

Plate Reader: Errors

🔥 **ReaderError(errorNumber, errorDescription)** - Error occurred during command execution.

| | | |
|------------------|---------|--------------------|
| errorNumber | Integer | Error number. |
| errorDescription | String | Error description. |

ReTiSoft Inc.
 366 Revus Avenue, Unit 21
 Mississauga, Ontario, Canada, L5G-4S5
 Main: 647-724-2398 Europe: 33-9-7518-0225
 Web: www.retisoft.ca Email: prodziew@retisoft.ca