



With ReTiSoft's Genera Suite, H-Cell and Staubli robotic arm integrated with a servo gripper you can have a flexible, robust and reliable automated system.

H-Cell complements our suite of software integration and scheduling products. It allows users to choose, reconfigure and reposition instruments on their automated systems.



Benefits

- High-end laboratory table with predrilled holes for instrument setup.
- Modular table with industrial grade components.
- Optional enclosure with temperature, humidity control and HEPA filters.
- TX40/TX60 Staubli robotic arm.
- SG0050 PTM servo gripper.
- Genera Suite integration framework and hybrid scheduler.
- Graphical representation of automated system workcell.
- Provides freedom from vendor-lock.

